ologenTM

Collagen Matrix for Tissue Repair

ologenTM Collagen Matrix
Operationstechnik

ologen™ Collagen Matrix- Operationstechnik

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Anhang I Referenzen

1. *Kohlhaas M.*, Trabeculectomy mit Mitomycin C und Ologen. Ophthalmol-Chirurgie 24: 19-24, 2012.
2. *Dada et al.* Trabeculectomy with combined use of subconjunctival collagen implant and low-dose mitomycin C. J. of Glaucoma 2012.
3. *Khoury A.S., et al.* Laser suture lysis through thick blebs using the Blumenthal lens. Arch.Ophthalmol. Vol.124, April 2006.
4. *Dietlein T.S., et al.* Secondary subconjunctival implantation of a biodegradable collagen-glycosaminoglycan matrix to treat ocular hypotony following trabeculectomy with mitomycin C. Brit. J. Ophthalmol, June 12,2013.
5. *Aptel F. et al.* Ultrasound biomicroscopy and optical coherence tomography imaging of filtering blebs after deep sclerectomy with a new collagen implant. Eur. J. Ophthalmol.19, 223-230, 2009.
6. *Shaarawy T. et al.* Biodegradable collagen implants in Trabeculectomy. Chapter 40, ISGS Textbook of Glaucoma Surgery, Jaypee, 2014.
7. *Rosentreter et al.*, Capsule excision and ologen™ implantation for revision after glaucoma drainage device surgery Graefes Arch Clin Exp Ophthalmol., 248:1319–1324, 2010.
8. *Oana S. and Vila J.* Tube exposure repair. J. of Current Glaucoma Practice, 6,3:139-142, 2012.
9. *Rosentreter A. et al.* Biodegradable implant for tissue repair after glaucoma drainage device surgery. Glaucoma, Vol. 23, no. 2, February 2012.
10. *Hamdi M.M and Hamdi I.M.* Scleral Repair by Biodegradable Collagen Implant in Strabismus Surgery. Ophthalmology Research: An International Journal, 3(4): 141-146, 2015, Article no.OR.2015.020.
11. *Cho E.Y. and Lee S.B.* Biodegradable collagen matrix (Ologen™) implant and conjunctival autograft for scleral necrosis after pterygium excision: two case reports. BMC Ophthalmology 15:140 DOI10.1186/s12886-015-0130-z. 2015.

Anhang II Studies of ologen™Collagen Matrix used in patients with 24 months follow-up or more.

I. *Cillino S. et al.* Biodegradable collagen matrix implant vs mitomycin-C as an adjuvant in trabeculectomy: a 24-month, randomized clinical trial. *Eye* (2011), 1–9.

II. *Cillino S. et al.* Biodegradable collagen matrix implant versus mitomycin-C in trabeculectomy: five-year follow-up. *BMC Ophthalmology* (2016) 16:24 DOI 10.1186/s12886-016-0198-0.

III. *El-Malah M.A. et al.* A prospective comparative study between trabeculectomy alone versus trabeculectomy with Ologen for patients with open angle glaucoma two years follows up. *Journal of American Science* 2014;10(11)

IV. *Fan K.C. et al.* Ologen versus Mitomycin-C for Trabeculectomy in a Predominantly African American Population. *Glauc Open Access* 2016, 1:1.

V. *Perez C.I. et al.* Trabeculectomy Combined With Collagen Matrix Implant (Ologen). *J. of Glaucoma*, DOI:10.109 JG.0000000000000551.2016.

VI. *Tanna A.P. et al.* Collagen matrix vs mitomycin-C in trabeculectomy and combined phacoemulsification and trabeculectomy: a randomized controlled trial. *BMC Ophthalmology* (2016) 16:217 DOI 10.1186 /s12886-016-0393-z.

VII. *Yuan F. et al.* Biodegradable 3D-Porous Collagen Matrix (Ologen) Compared with Mitomycin C for Treatment of Primary Open-Angle Glaucoma: Results at 5 Years. *Journal of Ophthalmology* Vol. 2015, Article ID 637537.

Anhang III Step-by-step description of a standard trabeculectomy with the ologen™Collagen Matrix, as performed in the Rotterdam Eye Hospital.

Accompanying text of the instruction video by Prof. Dr. Hans Lemij which may be found on <http://www.aeonastron.com/> and on <https://www.youtube.com/watch?v=Kybop3Zildo>.

- 0.01 A bridle suture is put through the cornea and the eye is gently pulled downward and fixed.
- 0.30 The cornea is covered with a sponge that will be moistened a little later to keep the cornea wet.
- 0.40 The conjunctiva and Tenon's capsule are opened with scissors and a limbus based flap is created.
- 0.59 Any adhesions are gently dissected.
- 1.17 Judicious cautery is applied to scleral vessels in the area where the scleral flap will soon be created. As little cautery as possible is performed to minimize the thermal damage that might possibly lead to undesirable scarring.
- 1.50 A scleral flap is created with a 30 degrees steel blade.
- 3.25 Additional cautery is applied to any leaking vessels.
- 3.44 A paracentesis is made, to fill the anterior chamber in later stages.
- 3.50 An entry into the anterior chamber is created under the scleral flap.
- 4.06 The trabeculectomy is performed by using a trephine.
- 4.20 A peripheral iridotomy is created, while an assistant holds the scleral flap.
- 4.35 Two 9-0 nylon sutures are put on the scleral flap, one on each corner. The bites into the surrounding sclera are made quite long to facilitate any future laser suturelysis. Sutures may be visualized with a Blumenthal lens. For instructions on how to use this lens see <http://www.glaucoma.org.il/physicians/volk-blumenthal-suturelysis-lens/intro>
- 5.20 The sutures are then tied and cut.

- 6.32 The anterior chamber is filled through the paracentesis to restore the intraocular pressure.
- 6.42 The knots of the sutures are then buried.
- 6.55 A 6 by 2 mm circular Ologen implant is then positioned onto the sclera, straddling the scleral flap and the neighboring sclera.
- 7.15 The Ologen implant soaks up aqueous that gradually effuses from the anterior chamber. The Ologen implant softens as it gets wet.
- 7.30 The conjunctiva and Tenon's capsule are then closed with a single running suture of 9-0 Vicryl.
- 9.42 At the end of the procedure, the anterior chamber is filled through the paracentesis and the bleb is gently formed and inflated.

Please visit www.ologen.com or
email info@aeonastron.com for any questions, comments
or order information.

Distributor Schweiz:

Sanassist AG

Hochwachtstrasse 34 | 6312 Steinhausen

www.sanassist.ch | info@sanassist.ch | 041 741 09 09

Aeon Astron Europe B.V. All rights reserved.

ologen® is a registered trademark owned by Body Organ Biomedical Corporation©